


MBT

Maxillary

Torque	-7°	-7°	0°	10°	17°	17°	10°	0°	-7°	-7°
ang.	0°	0°	8°	8°	4°	4°	8°	8°	0°	0°
width	2.8	2.8	2.8	2.8	2.8	2.8	2.8	2.8	2.8	2.8



Mandibular

Torque	-17°	-12°	0°	-6°	-6°	-6°	-6°	0°	-12°	-17°
ang.	0°	0°	3°	0°	0°	0°	0°	3°	0°	0°
width	2.8	2.8	2.8	2.8	2.8	2.8	2.8	2.8	2.8	2.8