

Elastomeric Products

Precision engineered and manufactured by G&H® in the USA



Actual patients, Julia and Samantha, wearing G&H elastomerics.

LIGATURE TIES/CHAIN



PEARL	GLOW BLUE	MET. PURPLE	FUCHSIA	YELLOW	GLOW	TURQUOISE	WHITE	
PEARL BLUE	MET. BLUE	LAVENDER	RUBY RED	MET. GOLD	PEARL GREEN	TRANSP. GREEN	MET. SILVER	BLACK
LT. SKY BLUE	DEEP SKY BLUE	PINK	MET. RED	ORANGE	NEON GREEN	TOOTH	SMOKE	
ROYAL BLUE	NAVY BLUE	FLUOR. PINK	MAROON	FLUOR. ORANGE	GREEN	CLEAR	GRAY	

*Dyna-Link™ Chain

Elastomeric Chain

Superior memory and durability with excellent color stability and stain resistance. Engineered to exacting standards to ensure consistent force delivery. **.026" Thickness** available in short, continuous and long links.

*Styles are not made with natural rubber latex.



Power-Linx™ Chain

A biomaterial breakthrough in elastic chain.

- Retains up to 85% of its original pull strength even after 6 weeks of engagement.
- Elasticity and recovery allows retention on any size bracket. **.020 Thickness** – Available in clear or gray only

*Bambino-Ties™ 1 **TOP SELLER!**

Soft Feel

12 Ties per Stick 84 Sticks/pkg (1,008 Ties total)

*Versa-Ties™ 2 **TOP SELLER!**

Firm Feel

22 Ties Per Stick 46 Sticks/pkg (1,012 Ties total)

*Select-Ties 3

Soft Feel

20 Ties Per Stick
50 Sticks/pkg (1,000 Ties total)

*Mini-Ties 4

Firm Feel

11 Ties Per Stick
92 Sticks/pkg (1,012 Ties total)

33
colors
to choose
from!



Elastomeric Ligature Color Display Booklet

Conveniently displays 33 elastomeric ligature tie colors for easy patient or clinician selection. Rotate color wheel to see how colors look together. Also shows all ligature tie part numbers for easy ordering/reordering. Measures 7½"x8½" when folded.

Private Logo Series

Printed with your logo!

Elastic carton contains
50 patient packs of 100 elastics.

Natural amber or
colors available.



Elastomeric Products



Precision engineered and manufactured by G&H in the U.S.A.

GHOrthodontics.com

TOLL-FREE 800-526-1026

orders@ghortho.com

G&H Orthodontics 2165 Earlywood Drive Franklin, IN 46131 U.S.A. 800-526-1026 U.S. & Canada +1 317-346-6655 International

MKT003-AM